

IMAGE CONTACT SHEET

High-resolution images available on USB drive.
Photos courtesy of Solvay.



1. Ultaire™ AKP Removable Partial Denture Frame

Image File Name: Ultaire AKP Frame

Removable partial denture frames made from Ultaire™ AKP offer a biocompatible, lightweight, metal-free, non-irritating and more aesthetically pleasing alternative to traditional metal RPD frames.



2. Ultaire™ AKP Removable Partial Denture Frame Without Teeth

Image File Name: Ultaire AKP Frame Without Teeth

Removable partial denture frames made from Ultaire™ AKP offer a biocompatible, lightweight, metal-free, non-irritating and more aesthetically pleasing alternative to traditional metal RPD frames.



3. Ultaire™ AKP Frame and Metal Frame

Image File Name: Ultaire AKP and Metal Frame

On the left: a removable partial denture frame made from Ultaire™ AKP, and on the right, a metal removable partial denture frame.

IMAGE CONTACT SHEET

High-resolution images available on USB drive.
Photos courtesy of Solvay.



4. Milled Dentivera™ Milling Disc

Image File Name: Milled Dentivera Milling Disc

A Dentivera™ milling disc made from Ultaire™ AKP illustrating three milled removable partial denture frames.



5. Dentivera™ Milling Disc Made from Ultaire™ AKP

Image File Name: Dentivera Milling Disc

The Dentivera™ milling disc made from Ultaire™ AKP is developed for milling removable partial dentures and is biocompatible, monomer- (MMA), BPA- and nickel-free. Ultaire™ AKP's base materials are currently used in permanent implants and non-implantable orthopedic, cardiovascular and renal devices.



6. Ultaire™ AKP Frame Held by Lab Technician

Image File Name: Frame Held by Lab Tech

This image shows a dental lab technician holding an Ultaire™ AKP removable partial denture frame.

IMAGE CONTACT SHEET

High-resolution images available on USB drive.
Photos courtesy of Solvay.

Workflow Comparison



7. Ultaire™ AKP Workflow Comparison

Image File Name: Workflow Comparison Graphic

The Ultaire™ AKP workflow increases lab efficiency by eliminating waxing, investing and casting procedures and reduces finishing and polishing time compared to metal.

Patient Advantages



8. Ultaire™ AKP Patient Advantages

Image File Name: Patient Advantages Graphic

Removable partial denture frames made from Ultaire™ AKP are tooth-supported, lightweight, biocompatible, leave no metal taste in the mouth and offer improved aesthetics versus metal.

9. Ultaire™ AKP Digital Workflow Time Savings

Image File Name: Time Savings Graphic

Dentivera™ milling discs made from Ultaire™ AKP are suitable for use in CAD/CAM digital workflows. Incorporating the digital workflow into business operations eliminates the waxing and casting steps of metal frames and provides approximately one hour in time savings per case, increasing productivity and ultimately profitability for dental labs.

